


---

# MEDIA RECIPES FOR PHENOTYPIC ASSAYS

**Created/updated by: Laura Rossi**    **Date: February 18, 2015**  
Surette Lab, McMaster University  
Hamilton, ON, Canada  
[www.surettelab.ca](http://www.surettelab.ca)



## BACKGROUND

- The following recipes allow for the characterization of bacteria with regard to its ability to produce the secreted proteins that are specified for each media type. References are given per media type where applicable.

## EQUIPMENT

- Scale, stir plate
- Autoclave

## BLOOD PLATES – assay the ability to break down blood cells – Hemolysis.

- Make CBA plates with 5% sheep blood.
- CBA agar is Columbia Blood Agar base, autoclave, cool and add 5% sheep blood.

## SKIM MILK PLATES – assay the ability to break down proteins – Proteinases.

- Make BHI agar with 1-1.5% final of skim milk. You can do this in any ratio you want but autoclave the two, BHI and skim milk, separately, then mix after autoclaving.
- Usually make 500mL of BHI agar and 20% skim milk. Autoclave, cool. Add 25mL of 20% skim milk to 500mL BHI for final of 1%.

## DNA PLATES – assay the ability to break down DNA – DNases.

- Make BHI agar, add 2g of DNA (Sigma#D1626) to 1L dH<sub>2</sub>O for a final of 2mg/mL DNA. Autoclave.
- After growth, for visualizing the zones, the plate needs to be flooded with 1N HCl, the zones will be visible within minutes.

### PRUSSIAN BLUE PLATES – assay the ability to produce $H_2O_2$ (Saito et al, 2007).

- Dissolve 1g  $FeCl_3 \cdot 6H_2O$  in 50mL  $H_2O$  and 1g potassium ferricyanide in separate 50mL  $H_2O$ . Mix the two solutions to give a reddish brown colour.
- Dissolve BHI agar for 1L into less than 900mL  $H_2O$ .
- Add the first mixed solution, 100mL, to BHI solution to form a green precipitate. Top up to 1L.
- Autoclave. Mix well before pouring plates.
- Zones will be bright blue where  $H_2O_2$  is being produced.

### HYALURONATE PLATES – assay the ability to produce hyaluronidases. (Rodney, Smith 1968)

- 3.7g BHI broth, 1g agar, 60mL  $dH_2O$ . Autoclave. Cool to 45°C-55°C, add 20mL of 5% Bovine Serum Albumin and 20mL of 2mg/mL umbilical cord hyaluronic acid. Pour plates. Cool. Refrigerate before use.
- For analysis flood plates with 2N acetic acid to visualize the zones.